

## ASTM D395 Method B Compression Set Fixture



Max Capacity: ~1,000 lbs.

Weight: 5kg each - steel

Spacer Dimensions: Includes 2x 6 Spacers with height of 4.5mm and 9.5mm.

**Notes:** This fixture is not used in conjunction with a UTM. The rubber sample is placed inside one of the slots and compressed to 25% of its original size. The spacers are used to prevent the fixture from compressing to more than 4.5mm or 9.5mm of sample height. Different spacer sizes are available. This fixture is designed to hold 4 of the larger standard specimens according to the ASTM drawing. Turn the screw to tighten the entire fixture and begin the long test.

## Schematic (mm)

